

## Software Release notes: Pippin Prep version 6.31 / 6.40

Current software version: Linux 6.31 / Windows 6.40

### Pippin Prep

Release notes:

Date: 08/07/18

Version: 6.31 / 6.40

### What's new from v.6.30 to 6.31 / 6.40

- Update works for both Windows and Linux based instruments
- 

Release notes:

Date: 6/1/17

Version: 6.31

### What's new from v.6.23 to 6.31?

- Progress Bar Update
  - Adjusted continuity test parameters
- 

Release notes:

Date: 9/1/16

Version: 6.23

### What's new from v.6.22 to 6.23?

Pippin Prep and Blue Pippin

- Software fix to eliminate occurrence of signal spike within the first few seconds of the run
- In peak-calling parameters, max allowable peak time increased to 50 mins

### Blue Pippin

- Updated 0.75% DF 50kb Marker Z1- Broad cassette definition
- Updated 0.75% DF 50kb Marker Z1- Tight cassette definition
- Added 3%DF Marker Q3 cassette definition - Size collection ranges 100-250bp
- Updated 0.75% DF 2-6kb Marker S1 cassette definition

### Pippin Prep

- Added 3%DF Marker P cassette definition - Size collection ranges 100-250bp
  - Adjusted continuity test parameters
-