

Software Release notes: Blue Pippin version 6.31 / 6.41

Current software version: Linux 6.31 / Windows 6.41

Blue Pippin

Release notes:

Date: 06/20/23

Version: 6.31 / 6.41

What's new from v.6.40 to 6.41

- Added New Range +T mode
-

Blue Pippin

Release notes:

Date: 08/07/18

Version: 6.31 / 6.40

What's new from v.6.30 to 6.31 / 6.40

- Update works for both Windows and Linux based instruments
 - Added New High Pass Plus Cassette Definitions
-

Release notes:

Date: 06/01/17

Version: 6.31

What's new from v.6.30 to 6.31?

- Progress Bar Update
 - Adjusted continuity test parameters
-

Release notes:

Date: 03/09/17

Version: 6.30

What's new from v.6.23 to 6.30?

This update is for BluePippin users only. Under heavy pulsed-field use with BluePippin, some older electrode arrays can form a salt bridge which may interfere with performance. This update runs a diagnostic routine to check for this. The routine is run automatically at the end of every run. The routine is not visualized by the user, but a warning screen will inform users if an error is detected.
